

Quantifit Use With PC or Stand-Alone Operation

Many of our customers asks us if it's better to use the Quantifit Respirator FitTester with a computer, or in "stand-alone" mode. The Quantifit will work well with either situation, but only you can decide which will work best for you.

As a rule of thumb: If you simply need to run tests, print them out, and be done with

it, then the stand-alone method is for you. If you are data-intensive and need to manage records year after year for a large number of employees, we would recommend operation through a computer. The in-between step would be to use the instrument in stand-alone mode and after the testing has been completed, transfer the data to FitTrack software.

	Operation Without Computer	Operation With Computer
Computer	No need for carrying around a laptop computer. Saves expense and risk of equipment loss.	Must have a computer connected to the Quanfitit.
Data	<p>Must use a designated HP printer to print calibration record and test results.</p> <p>Must input personnel information and test data for each subject.</p> <p>Can store up to 500 records onboard.</p> <p>Personnel information is limited to name and ID number.</p>	<p>Any Windows printer will work, and no need to print calibration record since it will be stored in the software database.</p> <p>If person has been previously tested, you can simply select their record and which mask they are using, and begin test.</p> <p>Can store an unlimited amount of records.</p> <p>Personnel information includes many options including custom information fields.</p>
Reports	<p>Reporting options include only a calibration report, fit test report, and fit test log.</p> <p>Generic report does not identify company.</p>	<p>Many reports available including sorting and filtering options, as well as the ability to export the data to Excel or a PDF file (pdf writer must be installed).</p> <p>Name of company will print at top of fit test report.</p>
Data Transfer	If you choose to operate the Quanfitit in stand-alone mode, you can transfer the data into FitTrack and receive all the benefits of FitTrack, and your records will merge with existing employee data.	

Öäcã c'áÁ[~]
 ÖÖÁQ c'!} æä } æÁQ&
 ÚUÁÓ[c'Á JI ÁÁO^T [c'ÁQ Á Î HFÉ
 Ì €€Í GÈGJHÁ/GFJÈÌ Ì È Ì G
 æÁCFJÈÌ Ì È Ì G
 • æ • O æ&ä dÈ { ÁÁ , , È&ä dÈ {